

# CLINICAL CHLOE - NURSING CARE SIMULATOR - OMNI2

PRODUCT NUMBER: 1260



## DESCRIPTION:

Clinical Chloe™ with OMNI® 2 is a full-body, adult-sized simulator designed for wireless training of general nursing clinical skills including tracheostomy procedures, bandaging, breast examinations, and more.

## SKILLS:

- Otic irrigation
- Patient mobilization with supine and sitting positions
- Bathing and bandaging
- Intramuscular injection exercises
- Male and female breast palpation exercises

## MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

## CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



Vaginal douching  
Pap smear test procedure  
Tracheostomy tube placement  
Nasal and oral intubation  
Blood pressure monitoring  
Gastric tube placement  
Stoma care  
Gastrostomy  
Enema administration  
Catheterization

## CASES AND PATHOLOGIES:

### OMNI 2 WIRELESS INTERFACE:

- The new OMNI® 2 features basic life support-specific simulation tools designed to help boost skill acquisition, enhance the simulated experience, and support assessment and debriefing. It's the one handy tool that can do it all.
- Wireless control at distances of up to 30 feet away
- Control over 35+ vital sign parameters including HR, ECG, RR, BP, SpO2, EtCO2, and much more
- Update vital sign values on-the-fly Trend vital sign improvements and deteriorations over time
- Comprehensive ECG library with 25+ preprogrammed rhythms
- Evaluate simulation participants individually or as a group
- Easily document provider actions with the built-in checklist
- Save, email, and print event logs with trend graphs for data-rich debriefing

## CHARACTERISTICS:

### GENERAL:

- Adult-sized, full-body patient simulator
- Specify light, medium, or dark skin tone at no additional charge
- Realistic eyes for ophthalmic exercises
- Simulated ear canal for otic drops and irrigation
- Upper and lower dentures for oral hygiene
- Soft, realistic face skin, hands, feet, fingers, and toes
- Articulating head, jaw, shoulder, elbows, wrists, knees, and ankles
- Wig for haircare exercises and surgical draping
- Supports bathing and bandaging activity
- Supports supine and sitting positions

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

#### **INJECTION TRAINING:**

- Intramuscular injection sites in deltoids, quadriceps, and upper gluteal

#### **BREAST PALPATION:**

- Interchangeable male and female chest inserts
- Breast palpation capability with 1 abnormal breast

#### **GYNECOLOGIC:**

- Vaginal douching and pap smear exercises
- Realistic vagina and cervix

#### **AIRWAY MANAGEMENT:**

- Tracheostomy tube placement
- Nasal and oral tube placement
- Gastric tube feeding and suction

#### **BLOOD PRESSURE AUSCULTATION:**

- Auscultate Phase I – Phase V Korotkoff sounds between systolic and diastolic pressures using real stethoscope
- Programmable auscultation gap
- Korotkoff sounds are automatically silenced if auscultation gap is selected
- Korotkoff sounds are responsive to the heart rate and the rate of cuff deflation
- Programmable blood pressure auscultation tutor

#### **BLOOD PRESSURE PALPATION:**

- Programmable, palpable radial pulse synchronized with heart rate
- Automatic return of radial pulse when cuff pressure is less than the selected systolic blood pressure

#### **GASTROINTESTINAL:**

- Transverse colostomy, ileostomy, and suprapubic stomas to practice irrigation
- Opening for gastrostomy
- Interchangeable male and female genitalia
- Realistic urethral passage and bladder for catheterization exercises
- Enema administration capability
- Detachable and removable internal tanks

## **INCLUDE:**

#### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

#### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- 1 Clinical Chloe manikin
- 1 OMNI 2 Wireless controller
- 1 Blood pressure training arm (installed)
- 1 BP cuff
- 2 Stomas
- 1 Power supply adapter
- 1 Carrying bag
- 1 Instructions manual

## ADDITIONAL MODULES:

### **ARTERIAL AND VENOUS TRAINING ARM (option):**

- A micropump delivers automatically generated arterial pulses at the radial and brachial sites and controls arterial blood flow by allowing variable heart rate and pulse strength.
- Interchangeable arterial and venous inserts within the forearm allow the creation of arteriovenous (AV) fistulas and for the placement of AV grafts, while a simulated healed fistula insert provides a platform on which hemodialysis exercises can be performed.
- An additional multi-layer insert in the bicep area can be used for incision and suture training exercises.

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>