

# EQUINE PALPATION AND COLIC SIMULATOR

PRODUCT NUMBER: 2483



## DESCRIPTION:

This simulator, modelled as a **standard 15 hand Quarter Horse** is perfect for teaching **palpation and exploration techniques** in veterinary students.

It also includes the equine **venipuncture and intramuscular** equine neck to extend training

## CHARACTERISTICS:

### HORSE BODY FEATURES:

- Modeled as a standard 15 hand Quarter Horse
- Steel reinforced epoxy/fibreglass construction with water resistant components for ease of cleaning
- 1.5m at the shoulder and 2.25m from nose to tail. Widest point is 0.6m

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

- 5 section inflatable latex GI Tract. Right and left ventral colon, right and left dorsal colon, and cecum. 3m section of small intestine with mount
- Spleen and left kidney attached with the renosplenic ligament. Palpable aorta capable of simulating a pulse
- Palpable soft uterus with cervix, urethra and broad ligament
- 3 sets of interchangeable ovaries; anestrus, estrus, and estrus with follicles
- Abdominocentesis function with replaceable pad; fluid reservoir and palpable linea
- Comes standard in black
- Custom markings available upon request
- Artificial, removable tail
- Soft, durable perineum panel
- Replica polyurethane pelvis
- Inflatable vinyl rectum
- Recycled plastic base

#### **EQUINE NECK VENIPUNCTURE AND INTRAMUSCULAR INJECTION AREA FEATURES:**

- Jugular venipuncture with palpable veins
- Easily replaceable latex vein tubing that can endure multiple punctures; includes 5 sets of veins
- Mechanical system allows the veins to be pressurized
- Gravity fed, no pump required
- Long lasting, easily removable, washable, replaceable hide covering with mane
- IM site that can endure multiple punctures and accept fluid injections
- Soft silicone head with flexible ears and well rendered anatomy illustrates facial injection sites and accepts a bridle
- Rolling Stand, adjustable for height, rostral tilt, and 90 degree rotation

#### **INCLUDE:**

- 1 Equine full-size simulator
- 3 Sets of interchangeable ovaries
- 1 Palpable soft uterus with cervix
- 1 Inflatable GI tract
- 1 Spleen and left kidney
- 1 Rolling stand
- 1 Removable leg (installed)
- 1 Instructions manual

#### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

#### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>