

Ultra-realistic adult male simulator - LifeCast

PRODUCT NUMBER: 3561



DESCRIPTION:

Ultra-realistic adult -male- simulator features realistic details, such as veins and underlying structures, hair, and an anatomically accurate airway for the correct training of various procedures.

They do not require electrical power to operate and are therefore ideal for on-site simulation sessions, which provide an effective and realistic extension to clinical practice-based training.

SKILLS:

- Realistic airway management

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Oral and nasal intubation
- Advanced Life Support (ALS)
- Cardiopulmonary resuscitation (CPR)
- Bilateral intraosseous infusion of tibia and humeral head
- Needle decompression in 2nd and 5th intercostal space
- Bladder catheterisation (*Partial)

CHARACTERISTICS:

GENERAL:

- Representation of an adult male, 40 years old and weighing 37 kg.
- Can be configured in a variety of ethnicities
- Created from real people using the life-casting and 3D scanning processes.
- Careful detailing to ensure immersion
- Mobility for manipulation
- No electrical connections required

INCLUDE:

- Ultra-realistic adult male simulator
- Instruction Manual
- Transport bag with wheels

ADDITIONAL MODULES:

- Integrated loudspeaker (Ref. 3582)
- Patient Monitor (Ref. 1010)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>