


Ultra-realistic female Senior simulator - LifeCast

PRODUCT NUMBER: 3581



DESCRIPTION:

Ultra-realistic -female- senior simulator, with realistic details, such as veins and underlying structures, hair and an anatomically accurate airway for the correct training of procedures. Created with a manageable weight and an average adult height.

They do not require electrical power to operate and are therefore ideal for on-site simulation sessions; they encourage more natural handling, bring a new level of realism to medical training, help create emotional attachment and commitment to clinical practice-based training.

SKILLS:

- Cardiopulmonary resuscitation (CPR)
- Staged airway management
- Advanced Life Support (ALS)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Bilateral intraosseous infusion in tibia and humeral head.
- Needle decompression in the 2nd and 5th intercostal space.
- Bladder catheterisation (*Partial)

CHARACTERISTICS:

GENERAL:

- Realistic airway anatomy
- Upper airway management

- Mouth, teeth and jaw
- Bag valve mask ventilation,
- Intubation (nasal and oral)
- Lower jaw manipulation
- Enables application of the chin manoeuvre to open the airway prior to ventilations
- Use of laryngoscopes and supraglottic devices among others.

- Compatible with CPR

- Its internal structure reproduces the reality of a CPR manoeuvre on an elderly person.

- Nostrils for nasopharyngeal airway
- Intraosseous insertions in the humeral head and bilateral tibia.
- Decompression needle insertion: 2nd intercostal space (frontal) and 5th intercostal space.
- Standard weight 25 & 30 kg (estimated)**.
- Standard age approximately 80 years
- Closed eyelids capable of manual opening
- Perforated hair and body hair (Genitalia)

INCLUDE:

- Ultra-realistic female geriatric simulator
- Manual de Instrucciones
- Bolsa de transporte con ruedas

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



ADDITIONAL MODULES:

- Integrated loudspeaker (Ref. 3582)
- Patient Monitor (Ref. 1010)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>