

Ultra-realistic 3 year old child simulator **LifeCast**

PRODUCT NUMBER: 3607



DESCRIPTION:

Ultra-realistic 3 year old child simulator, with realistic details, such as veins and underlying structures, hair and an anatomically accurate airway for the correct training of procedures. Created to reflect the challenges posed by children of this age, and the procedures associated with emergencies.

Created with a manageable weight and average adult height. They do not require electrical power to operate and are therefore ideal for on-site simulation sessions; bringing a new level of realism to medical training, they help to create emotional attachment and commitment to clinical practice-based training.

SKILLS:

- Realistic airway management

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Oral and nasal intubation
- Cardiopulmonary resuscitation (CPR)
- Advanced Life Support (ALS)

CHARACTERISTICS:

GENERAL:

- Representation of an 3-year-old boy
- Measurements: 80 cm and 14kg
- Can be configured in a variety of ethnicities
- Created from real people through the processes of life-casting and 3D scanning
- Careful detailing to ensure immersion
- Mobility for handling
- No electrical connections required

INCLUDE:

- Ultra-realistic 3 year old boy simulator
- Instruction Manual
- Transport bag

ADDITIONAL MODULES:

- Bilateral intraosseous infusion in tibia (Ref.3608)
- Patient Monitor (Ref. 1010)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>