


# MIKE MICHELLE - NURSING CARE PEDIATRIC 5 YEAR - OMNI2

PRODUCT NUMBER: 1274



## DESCRIPTION:

The **Five-Year-Old PEDI® Nursing Care Patient Simulator** is designed to help participants practice and develop pediatric nursing care and resuscitation skills.

It includes an intubatable airway, an injection training arm, an intraosseous injection leg, and appropriate arterial/venous system to ensure participants in ALS training programs receive hands-on training.

## SKILLS:

- Oral and nasal endotracheal intubation

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- NG feeding
- Intraosseous infusion
- Intramuscular injection exercises
- IV training
- Manual palpation of pulses
- BVM ventilation
- Resuscitation training (CPR)
- Lavage/gavage
- Ear canal examination
- Gastric suctioning techniques

## CHARACTERISTICS:

### FEATURES:

- Full-body simulator with soft, lifelike face skin and molded hair
- Eyes can be closed or opened manually
- Head tilt/chin lift
- Endotracheal intubation
- Nasogastric intubation
- Tongue, teeth, and realistic airway
- Intraosseous infusion
- IV training arm
- Femoral vein
- Manual palpable pulses: carotid, femoral, brachial
- BVM with realistic chest rise
- CPR anatomical landmarks
- Lavage/gavage
- Intramuscular injection sites
- Simulated ear canal
- Practice gastric suctioning exercises
- Articulating head, jaw, arm, and legs

## INCLUDE:

- 1 Multipurpose Pediatric patient simulator
- 1 OMNI2 Tablet controller
- 1 Squeeze bulb with tube
- 1 IO leg tibia bones

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



**More Than Simulators™**  
*Excellence. Every patient. Every time.*

- 1 Right tibia bone cover
- 1 Right arm skin
- 1 Filling syringe
- 1 Synthetic lubricant
- 1 Blood concentrate
- 1 Instructions manual



**MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

**CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>