


FACIAL INJECTION SIMULATOR - MUSCLES, ARTERIES AND NERVES

PRODUCT NUMBER: 3180



DESCRIPTION:

This facial injection simulator is made of a high-quality platinum-hardened prosthesis silicone with a realistic tactile feel. The skins are then mounted on a rigid foam skull that has various landmarks to aid the teaching process. These skulls are attached to 2 flexible stands that allow the head to be tilted into suitable positions. Each head is then hand painted, attached to a plastic backing, and mounted on a high quality base.

This version shows almost all muscle groups, the important arteries and nerves on one side. It is perfect for training because you can explain the anatomy and show you where to inject correctly.

SKILLS:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Botox (botulinum toxin) injection
- Superficial filling techniques
- Facial anatomy recognition

CHARACTERISTICS:

GENERAL:

- Made from ultra-realistic high quality silicone
- Expertly sculpted and moulded
- Hand-painted down to the smallest detail

INSTRUCTIONS FOR USE:

- Allows fluid injection: After administration, dry the treated area with a tissue and/or paper. Pinch the area to remove the fluid from the inside through the same hole where the injection was made.
- Check that the fluid has been completely removed to avoid deformation of the appearance of the simulator.

INCLUDE:

- 1 Facial injections model
- 1 Stand

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>