


# FACIAL INJECTION SIMULATOR - ADULT F-MODEL

PRODUCT NUMBER: 3183



## DESCRIPTION:

Our facial injection simulator is made of a high quality platinum-hardened silicone prosthesis with a realistic tactile feel. The skins are then mounted on a rigid foam skull that has several reference points to aid in the teaching process.

These skulls are attached to 2 flexible supports that allow the head to be tilted into suitable positions. Each head is then hand painted, attached to a plastic backing and mounted on a high quality base.

## SKILLS:

- Botox injection (botulinum toxin)
- Superficial filler techniques

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



## CHARACTERISTICS:

### GENERAL:

- Made from ultra-realistic high quality silicone.
- Expertly sculpted and moulded
- Hand painted down to the smallest detail

### INSTRUCTIONS FOR USE:

- Allows fluid injection: After administration, dry the treated area with a tissue and/or paper. Pinch the area to remove the fluid from the inside through the same hole where the injection was made.
- Check that the fluid has been completely removed to avoid deformation of the appearance of the simulator.

## INCLUDE:

- 1 Model of facial injections
- 1 Support

## ADDITIONAL MODULES:

### MODELS AVAILABLE:

- Model A (Ref. 3182): young woman with fair skin.
- Model B (Ref. 3179): Older woman with fair skin.
- Model C (Ref. 3181): Dark-skinned middle-aged woman
- Model D (Ref. 3178): Elderly man with fair skin
- Model E (Ref. 3177): Young light-skinned woman
- Model F (Ref. 3183): Middle-aged woman with medium skin
- Model G (Ref. 3175): middle-aged woman with fair skin

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>