



Adult intubation head - AirSim

PRODUCT NUMBER: 3443



DESCRIPTION:

The adult intubation head trainer offers a solution for varying degrees of difficulty when practising airway management.

This versatile model ranges from a routine intubation scenario to the most difficult and complex intubation scenario; facilitating a variety of training scenarios for different levels of practitioners, from beginner to advanced.

The advanced airway management training head is ideal for training in intubation, percutaneous tracheostomy and surgical and needle cricothyroidotomy.

SKILLS:

- Laryngospasm

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Tongue oedema
- Airway management
- Nasal, oral and digital intubation
- Needle and surgical cricothyroidotomy
- Percutaneous tracheostomy
- Double nasotracheal intubation
- Bag and mask ventilation techniques
- Full range of supraglottic devices
- Direct laryngoscopy
- Endotracheal tube insertion
- Awake fiberoptic examination
- Combined tube insertion
- Videolaryngoscopy (including Sellick manoeuvre)
- Nasogastric tube insertion techniques

CHARACTERISTICS:

GENERAL:

- Laryngospasm
 - Adjustable from normal to complete airway obstruction.
 - Sudden vocal cord closure can be simulated.
 - When successfully acted upon by the user, the laryngospasm can be mechanically reversed to provide positive feedback.
- Displaced larynx
 - Ability to move the larynx to a more anterior position to make laryngoscopy more difficult.
 - Adjustable to provide various levels of difficulty.
- Tongue oedema
 - The tongue can be swollen to simulate various degrees of macroglossia.
- Retracted/protruding jaw
 - Ability to simulate 2 degrees of mandibular protrusion; a normal jaw and a retracted jaw.
 - Simulated overbite and underbite characteristics

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Trismus
- Ability to coarsely restrict mouth opening to make insertion of airway devices difficult
- Neck locking/head rotation restriction
- Possibility to simulate cervical arthropathy by restricting range of motion of the cervical spine
- Neck block/head rotation restriction
- Airway management at the front of the neck
- Simulating needle and surgical cricothyroidotomy and percutaneous tracheostomy
- Realistic, robust, innovative and high quality intubation simulator.

ANATOMÍA COMPLETA:

- Tooth
- Tongue
- Oral and nasal pharynx
- Larynx
- Epiglottis
- Arytenoids
- False vocal cords
- True vocal cords
- Trachea

INCLUDE:

- 1 Adult intubation head
- 1 carrying case
- 100 ml lubrication bottle
- 5 larynx inserts
- 1 neck skin wrap
- 2 sets of adhesive pads (attachable to suction cups to increase stability if needed)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>