


SkinPad - Incision and suture skills trainer

PRODUCT NUMBER: 1017



DESCRIPTION:

Completely new **3 layer skin pad** for demonstrating and practising a wide variety of incisions and suturing skins.

It is available in **medium size** (120x75mm) or **large size** (120x150mm).

SKILLS:

- Incisions: linear, ellipse, flaps and shaped
- Intradermal suturing
- Simple and advanced interrupted suturing

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Subcuticular suturing
- Continuous suturing
- Stapling
- Use of adhesive strips

CHARACTERISTICS:

- Dimensions: 120 x 75 mm o 120 x 150 mm
 - Realistic representation of the 3 main skin layers: epidermis, dermis and subdermal layer (fascia).
 - Realistic retention of sutures thanks to its high-quality silicone
 - Integrated mesh to increase the resistance of the tissue
 - Adequate for beginners and specialists
 - Reusable. The entire surface allows incisions and anaesthetic injection.
- It is recommended to use for one year after purchase. High quality silicone includes oils in its manufacture that guarantee a correct texture. These oils evaporate over time.

INCLUDE:

- 1 SkinPad Thin (can be ordered in packs of 2, 4, 6 or 10 units)

ADDITIONAL MODULES:

- **INTRUMENTS CASE (option):** Add a 19-pieces dissection kit to enhance your suture practice.
- **LIPOMA PAD (sold separately):** A realistic pad with 2 lipoma to practice full lipoma removal.
- **SEBACEOUS CYST PAD (sold separately):** A realistic pad with 3 cysts to practice full sebaceous cyst removal procedure.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>