

Adult CPR simulator with indicator light - PRESTAN

PRODUCT NUMBER: 2967



DESCRIPTION:

The PRESTAN adult simulator is realistic to the eye and touch, like no other CPR training simulator on the market. Lightweight and quick and easy to set up. Equipped with a revolutionary CPR rate monitor that provides instantaneous rate and depth information for both instructor and trainee.

CHARACTERISTICS:

GENERAL:

- Includes the latest AHA Integrated Feedback Directive.
- Available in medium or dark skin tones, CPR feedback only and in single, 4-unit or 4-unit diversity kit configurations.
- All simulators are available with a revolutionary CPR rate monitor that provides instant feedback for both instructors and trainees. The visual feedback from the CPR rate monitor gives trainees a true sense of performing 100-120 compressions per minute. Lights continuously signal rate progression, confirming that compressions are within the 100-120 cpm range recommended by industry guidelines. The yellow light will flash if the upper limit of 120 cpm is exceeded, warning trainees to slow down. Instructors can monitor multiple trainees quickly and easily.
- In addition to the visual CPR rate monitor, trainees will also hear a clicking sound when the chest is pushed to the proper depth of 2.0 to 2.4 inches, allowing them to experience the true force required to perform actual chest compressions.

INCLUDE:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



More Than Simulators™
Excellence. Every patient. Every time.

- Prestan 2000 Manikin
- Hygienic bags for mouth-to-mouth ventilation
- User manual
- Carrying bag



MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>