

CPR SENSOR WITH FEEDBACK - SKILLQUBE CPR

PRODUCT NUMBER: 3261



DESCRIPTION:

The sensors provide the instructor with information about the pressure depth and speed in the integrated display on the scenario controller.

The arrows above and below the respective values indicate whether, for example, the button is pressed too quickly or too deeply.

CHARACTERISTICS:

GENERAL:

- Thanks to the intelligent sensor recognition, even with larger teams, only the sensor that is also performing cardiopulmonary resuscitation is displayed in the controller.
- In the controller, the instructor gets to a more detailed menu in which the last two minutes are graphically displayed in order to be able to better assess the progress and quality of cardiopulmonary resuscitation.

READINGS:

- Ratio of successful and unsuccessful compressions in%
- CPR fraction in%
- Average team performance
- The respective individual evaluations can be viewed in the logbook

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



More Than Simulators™
Excellence. Every patient. Every time.

INCLUDE:

- 1 CPR intelligent sensor
- 1 Software license (1 year)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>