



# Nasogastric Probe Trainer

PRODUCT NUMBER: 1469



## DESCRIPTION:

This simulator allows to learn how to place the NG tube. The oesophagus and stomach are stored in the head for easy storage and transport from classroom to classroom or across laboratory tables. The head is mounted on a base that allows it to stand securely. The absence of moving parts and simple construction ensure longevity of the product. The only interchangeable part is the stomach pouch. Cleaning is quick and easy.

## SKILLS:

- Nasogastric tube care and insertion training
- Stomach elevation

## CHARACTERISTICS:

### GENERAL:

- Weight: 1.35kg
- Dimensions: 30.5x30.5x23 cm

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>