

# HAL - SURGICAL AIRWAY & CPR SKILLS TRAINER

PRODUCT NUMBER: 1222



## DESCRIPTION:

HAL® is a portable airway management and CPR skills trainer. HAL was designed to help you teach basic and advanced airway management as well as cardiopulmonary resuscitation. HAL helps participants learn to recognize laryngospasm and practice performing emergency tracheotomies and cricothyrotomies.

HAL® is available in a full-body configuration.

## SKILLS:

- Airway management
- Nasal and oral intubation
- Resuscitation training (CPR)

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- BVM ventilation
- Emergency airway management (tracheotomy or cricothyrotomy)
- Detection of laryngospasm

## CHARACTERISTICS:

### FEATURES:

- Adult-sized, male upper torso
- Anatomically correct landmarks for proper CPR hand placement
- Manually open or close eyes
- Realistic head tilt/chin lift and jaw thrust
- Durable, lifelike airway anatomy
- Nasal and oral Intubation: ETT, LMA, King LT
- Laryngospasms
- Supports tracheostomies and cricothyrotomies
- Palpable carotid pulse
- Visible chest rise with bag-valve-mask ventilation
- Gastric distension with esophageal intubation
- Right mainstem intubation presents unilateral chest rise
- Realistic chest compression recoil and depth

## INCLUDE:

- 1 HAL Skills trainer
- 1 Carrying case
- 1 Instructions manual

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>