

# CHEST TUBE MANIKIN

PRODUCT NUMBER: 2622



## DESCRIPTION:

This manikin is designed specifically to **teach the theory, anatomy, and skills needed to manage pre-hospital chest trauma**, as well as ongoing chest tube maintenance.

## SKILLS:

- Pre-hospital chest trauma complete management
- Chest tube surgical placement and maintenance for pleura effusion
- Tension pneumothorax drainage
- Human chest anatomical structures teaching

## CHARACTERISTICS:

- The right side of the manikin has two cut-away viewing areas to provide awareness of the anatomical relationships between the skin surface, musculature, ribs, and lungs.
- The left side has a pressurized tension pneumothorax site to relieve air that has accumulated within the pleural space and is restricting lung inflation.

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

- There is also a site where chest tubes may be surgically placed to treat pleural effusion by draining fluids from the pleural space.
- The fluid color, volume, and viscosity is controlled by the instructor.

## INCLUDE:

- Chest tube manikin
- 5 Replaceable surgical chest tube site pads
- 5 Replaceable pneumothorax chest pads
- 1 Manual pump
- 1 Fluid bag reservoir
- 1 Hard carry case
- Instructions manual

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>