

Small vein simulator

PRODUCT NUMBER: 3465



DESCRIPTION:

This two-layer vein offers the opportunity to train surgical ligatures as well as different types of anastomosis.

This two-layer vein simulator is especially thin, perfectly imitating the characteristics of the delicate venous tissue of dogs and cats. It is suitable for practising incisions, sutures and different anastomosis techniques under realistic conditions. The visual and textural demarcation between the tunica intima and the tunica externa offers the opportunity to practice various sutures.

SKILLS:

- Incisions
- Sutures
- Anastomotic techniques

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

CHARACTERISTICS:

GENERAL:

- Reusable several times
- Tear-resistant and retain their shape and texture
- Practical handling with the vein holder for vertical or horizontal positioning as required.
- Length: 24mm
- Diameter: 6mm
- Layers
- Inner tunic
- Outer tunic

INCLUDE:

- 1 Small vein simulator

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>