




# BUCKLE - Retinal implant

PRODUCT NUMBER: 1817



## DESCRIPTION:

Practice **buckling techniques** with realistic sclera and rectus muscles.

Hollow model allows simulation of **retinal implant techniques** in an open-sky platform

## SKILLS:

- Buckling techniques
- Retinal implan techniques
- Instrument manouvering skills

## CHARACTERISTICS:

### KEY FEATURES:

Realistic sclera and rectus muscles

Open-sky advanced platform

Hollow model to allow simulation of retinal implant techniques.

### INSTRUCTIONS FOR CARE:

## MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

## CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- \* Store in a cool, dry and dark place (a drawer will be fine). Extended exposure to sunlight (UV) may affect material properties. Prolonged exposure to humidity or high temperatures may adversely affect material properties.
- \* Do not place heavy objects on top of the model. Prolonged compression may deform the model. Store with the lens down in its hardbox.
- \* Clean with anti-bacterial soap and water. Do not soak the models for cleaning. Avoid using alcohol or bleach. Gently scrub with cotton swab if necessary.
- \* Air or blow dry before storing. Do not store wet.

## INCLUDE:

- 2 Hollow globe eye simulator
- 1 Storage box

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>