

OSCAR - DIDACTIC SKELETON

PRODUCT NUMBER: 1998



DESCRIPTION:

First class reproduction of a male adult skeleton. This skeleton is manufactured after our strict quality regulations.

The stable mounting and robust structure guarantees a long lasting quality even in intensive use.

CASES AND PATHOLOGIES:

Herniated disk

CHARACTERISTICS:

- Made of durable, unbreakable plastic

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Natural skeleton size
- The teeth are manufactured separately and inserted into the jaw.
- The limbs can be removed and reattached very fast and easy using a quick fastener.
- Shoulder, hip and knee joints do have flexible rubber attachments that allow lifelike, sliding joint movements. This distinguishes Oscar pleasantly from less mobile models with metal joints.
- The about 200 bones of the skeleton correspond with real human bones in size and nearly in weight.
- The skeleton is anatomically correct and complete, it shows all important structures and foramina.
- 5-spoke roller stand, rugged, smooth-running and removable.
- The spine of the skeleton shows spinal nerves, the vertebral artery and a dorsolateral disc prolapse.
- 3 year warranty, lifelong repair service
- Size with stand: 178cm
- Size without stand: 160cm
- Weight: 9 kg

INCLUDE:

- 1 Didactic skeleton
- 1 Stand

ADDITIONAL MODULES:

DUST COVER FOR SKELETON MODELS (option):

- This dust cover is made of sturdy plastic fabric, it is tear resistant and provides perfect protection against dust accumulation.
- The cover is simply pulled over the skeleton and is virtually indestructible because it has no fragile parts like e.g. zippers that could be destroyed.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>