


BOVINE VASCULAR ACCESS SIMULATOR

PRODUCT NUMBER: 2769



DESCRIPTION:

The **Bovine Injection Simulator** is a one-of-a-kind education tool. Use this simulator to **teach proper cattle injection techniques** and ear tagging as well as growth implant basics.

SKILLS:

- Cervical injections
- Intravenous injections
- Intramuscular injections
- Subcutaneous injections
- Piercing placement
- Implant placement

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



CHARACTERISTICS:

GENERAL:

- A full-size trainer for cattle injection is a unique teaching tool with versatile and multitasking application.
- The vascular access simulator is engineered to provide the most anatomically correct and realistic simulated venipuncture training experience available.
- The vascular access simulator delivers excellence in veterinary training.
- Dimensions: 71x102x173cm
- Weight: Approx. 18kg

INCLUDE:

- Cow head with neck
- Removable training cover
- Drain system for liquids
- Stand with wheels
- Foot pump for pulse imitation
- 10 Pairs of exchangeable ears

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>