


# Bovine coccygeal vein simulator

PRODUCT NUMBER: 3795



## DESCRIPTION:

Model of a cow's tail with artificial blood pump and realistic coccygeal vein for practising blood collection techniques.

## CHARACTERISTICS:

### GENERAL:

- Anatomically correct
- The vertebrae and the position of the veins have been reproduced from real bovine anatomy. real bovine anatomy, allowing students to palpate the correct.
- Blood pump
- The blood is circulated through the vein by a small electronic pump which simulates the pumping of the heart.
- Needle resistance
- Realistic resistance is felt under the skin when the needle is inserted.

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

- Bulging vein
- When pressure is applied, the vein bulges so students can see where to insert the needle.
- Universal holder
- Supplied with a sturdy metal stand that attaches to a bench, allowing the tail to hang down in a natural position.

## INCLUDE:

- 1 Blood sampling tail with insert for coccygeal vein
- 1 Blood pump
- 1 Bottle of artificial blood

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>