


Equine Neck for venipuncture and IM injection

PRODUCT NUMBER: 2484



DESCRIPTION:

The **Equine Neck Venipuncture and IM injection simulator** is a stand-alone task trainer molded from a real horse that allows students to develop their veterinary skills on equines.

CHARACTERISTICS:

GENERAL:

- Easily replaceable latex vein tubing that can endure multiple punctures
- Mechanical system allows the veins to be pressurized
- Long lasting, easily removable, washable and replaceable mane
- Soft silicone head with flexible ears and well rendered anatomy illustrates facial injection sites and accepts a bridle
- Rolling stand with adjustable height, rostral tilt and 90 degree rotation

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



VENIPUNCTURE AND INJECTION:

- Jugular venipuncture with palpable veins
- IM site can accept fluid injections

INCLUDE:

- Equine venipuncture neck simulator
- 5 Sets of veins
- 1 Rolling stand

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>