


Shoulder simulator for ultrasound-guided infiltration

PRODUCT NUMBER: 3973



DESCRIPTION:

Shoulder simulator for infiltration suitable for ultrasound, representing a full-scale human shoulder.

SKILLS:

- Ultrasound-guided access
- Infiltration treatment

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

CHARACTERISTICS:

- The shoulder simulator contains the following structures:
 - Humerus
 - Scapula
 - Clavicle
 - Acromioclavicular ligament
 - Biceps ligament.
- The shoulder can be punctured with 27 (0.4 mm) to 22 (0.7 mm) gauge needles.

INCLUDE:

- 1 Shoulder simulator for ultrasound-guided infiltration
- 1 Base

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>