

THREE-VEIN VENOUS ACCESS TRAINER

PRODUCT NUMBER: 3166



DESCRIPTION:

This three-vein venous access trainer easily attaches to a person's arm so that the student can practice venous puncture techniques on a live subject.

It allows proper positioning on a real human arm, without the risks associated with puncturing live tissue. The Venatech IV trainer is ideal for practicing the initiation of IVs in a moving ambulance.

CHARACTERISTICS:

GENERAL:

- The trainer includes the three main veins used in phlebotomy: cephalic, basilic and median ulnar, all in the correct anatomical position.
- An arterial vessel is present in the trainer for simulation of arterial blood gases.
- It can also be used as a tool to show how to correct puncture through a vein into an artery.
- The trainer helps students learn essential communication and patient interaction skills.
- Dimensions: 17 x 8 x 3 cm
- 5 year warranty

INCLUDE:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



More Than Simulators™
Excellence. Every patient. Every time.

- 1 Venous access trainer
- 1 Arm band
- 2 Syringes for fluid insertion



MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>