



- NG exercises
- Umbilical catheterization
- Intraosseous infusion
- Injection and intravenous exercises
- Palpation of pulses
- Lavage/gavage
- Debriefing exercises

CARACTERISTICAS:

FEATURES:

- Full-body, full-term infant with articulating head, arms, legs
- Realistic airway with tongue, vocal cords, trachea, and esophagus for airway management exercises
- Support oral or nasal intubation plus suctioning using standard devices
- Realistic lung expansion with bag valve mask ventilation
- Programmable peripheral and central cyanosis as well as healthy skin tone
- Intravenous (IV) access on right arm and lower left leg support infusion and bolus injections
- Right tibia intraosseous access and infusion
- Palpable pulses: right brachial, radial, femoral and umbilical (manual bulb)
- Practice umbilical catheterization and infusion
- Simulate neonatal CPR with either two-thumb "encircling" technique or two finger alternate compression method
- Supports gastrostomy tube placement

INCLUI:

- 1 PEDI Blue neonate simulator
- 1 OMNI2 Wireless interface (tablet+license)
- 2 Replacements for tibia IO access
- Power supply adapter
- 1 Carrying case
- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>