

TRUECLOT - Simulated Hemostatic Gauze Z-fold (3 unit)

Nº DE PRODUTO: 2652



DESCRIÇÃO:

Simulated Haemostatic Gauze is designed to work in conjunction with our TrueClot® Blood Simulant to create realistic simulated blood clots.

CARACTERÍSTICAS:

- Each gauze package contains one strip of gauze measuring 3 inches wide and 12 feet long and is packaged Z-folded.
- They are packaged in a heat-sealed, white plastic pouch that is very similar to QuikClot® gauze packaging.
- TrueClot® Simulated Hemostatic Gauze is for simulation training only. It is not sterile and is not intended for medical use.

INCLUI:

- 3 Training Hemostatic gauzes, individually packed.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>