


BUDDY PLUS PACK II - ADULT/INFANT CPR TRAINING PACK WITH APP

Nº DE PRODUTO: 2941



DESCRIÇÃO:

Dramatic, realistic, affordable, and **2019 AHA Compliant** simulation package for **teaching BLS an CPR skills**.

CARACTERÍSTICAS:

GENERAL:

- Enables accurate and effective training with apps and sensor kit
- Student and instructor apps provide real-time feedback and feature adult capabilities, as well as a realistic interface
- **Basic Buddy®** manikin simple to assemble and maintain
- Provides an open airway using the head-tilt/chin-lift method with a visible chest rise when ventilated

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>



- Offers an easy-to-simulate obstructed airway
- Xiphoid process provides anatomical reference point for hand placement and compressions
- Training data is saved digitally on the phone or tablet
- Use only AED training devices. When using foam-bodied manikins for AED training, use only pads designed for foam bodies.

HEARTISENSE STUDENT AND INSTRUCTOR APP:

- Student and instructor apps run on Apple® iPhone® devices.
- Apple® iPad® version apps will be available February 2019.
- Android™ version apps will be available in the first half of 2019. (iPad®, Android™, and Apple® tablet/phone not included).

INCLUI:

- 4 **Basic Buddy®** adult/child manikins with head
- 40 **Basic Buddy®** lung/mouth protection bags
- 4 Adult foam compression blocks (for use with manikin only)
- 6 Lung insertion tools
- 2 **Baby Buddy®** infant manikins
- 20 **Baby Buddy®** lung/mouth protection bags
- CPR Add-on kit student* and instructor* apps
- 4 Add-on kits
- 4 Breath modules (each breath module includes lung/mouth protection bags with breath module holders)
- Instruction manuals

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>