

ADULT TORSO FOR ULTRASOUND TRAINING

Nº DE PRODUTO: 2748



DESCRIÇÃO:

Adult Human Torso (See-through/Transparent Version)™ is a realistic and anatomically correct phantom. It has been developed to provide simulated training for multiple Sonography or ultrasound courses.

This phantom can be customized with any pathology for learning the techniques of **FAST** (focused assessment with sonography for trauma), **FAMUS** (Focused Acute Medicine Ultrasound) and any other procedural training (Pericardial effusion, etc...)

CARACTERÍSTICAS:

FEATURES:

- It is a replica of a normal adult human being. An ideal phantom for training ultrasound technicians and other medical students.
- The model has the following anatomical structures:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>



- Heart
- Lungs
- Liver
- Stomach
- Diaphragm
- Spleen
- Pancreas
- Kidneys
- Vena cava
- Esophagus tube
- Aorta
- Colon
- Bladder
- Large intestine
- Small intestine
- Pelvic
- Ribcage
- Spine

GENERAL:

- Soft tissue and organs: Composition of urethane base soft resin
- Synthetic bones: Urethane base hard resin
- Size: 85 x 46 x 26 cm(approx.)
- Weight: 40 Kgs (approx.).

INCLUI:

- 1 Ultrasound-compatible adult torso trainer
- 1 Design video
- 1 Instructions manual
- 1 Box

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>