


ADV ZACK - Advanced multipurpose male care

Nº DE PRODUTO: 2910



DESCRIÇÃO:

The ZACK™ Multipurpose Male Care Simulator is a **full-sized, adult male lower torso** (abdomen and pelvis) designed to provide health care professionals with hands-on **training in the skills required to perform men's healthcare procedures**.

It is useful for demonstrating and practicing catheterization exercises, prostate examinations, testicular examinations, no-scalpel vasectomy, and rectal/colon examinations with life-like realism.

HABILIDADES:

- Catheterization exercises
- Prostate examination
- Testicular examination

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>



- Vasectomy
- Colonoscopy
- Rectum examination

CASOS E PATOLOGIAS:

COLON PATHOLOGIES:

- One melanoma
- Two ulcerative adenocarcinomas
- Adenomatous polyp
- Two sessile adenocarcinomas
- One advanced polypoid adenoma
- One pedunculated adenomatous polyp
- One pedunculated polyp

PROSTATE PATHOLOGIES:

- Enlarged benign mass
- Two discreet nodules
- Large mass
- Malignant invasive cancer

CARACTERISTICAS:

GENERAL:

- Adult male lower torso with lifelike feel and appearance
- Three available skin tones
- Easily replaceable vas deferens tubing allows for the repeated practice of vasectomy procedures
- Male genitalia with a patent urethra and rectum
- Removable lid for easy assembly and disassembly and access to the internal hardware
- Urinary bladder

PROSTATE EXAMINATION:

- Learn vital rectal examination and palpation skills needed to distinguish between the "feel" of a normal prostate, benign conditions, and

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>



various grades of prostatic malignancies.

- Four prostate inserts to help develop haptic skills required during intimate examinations

DIAGNOSE VARIOUS PROSTATE PATHOLOGIES:

- The internal push-button mechanism allows you to rotate through and position four different prostate inserts (enlarged benign mass, two discreet nodules, a large mass, and malignant invasive cancer) during digital rectal examination training.

COLONOSCOPY:

- Use a real endoscope to identify benign and malignant masses in the rectum and colon; practice biopsy and polyp removal with non-electrocautery devices.

TESTICULAR EXAMINATION:

- Realistic testes allow participants to develop palpation skills needed for intimate examination procedures.
- Four easily detachable testes (two normal and two abnormal) with distinct masses that can be identified during palpation

INCLUI:

- 1 Lower torso male simulator
- Scrotum
- Normal testes
- Rotational prostate assembly
- Cancerous testis 1
- Cancerous testis 2
- Fill and drain Kit
- Talcum powder
- Mineral oil
- 1 Soft carrying bag
- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>