



FULL TERM BIRTHING BABY ROHAN

Nº DE PRODUTO: 2053



DESCRIÇÃO:

Baby Rohan is a full Soft silicone Sim baby with a super realistic newborn look.

HABILIDADES:

- Birthing and delivery exercises
- Normal or emergency C-section training
- Vaginal birthing in many models of birthing manikins
- Intubation and ventilation
- Blood draw from umbilical vein
- New parenthood training; bathing and dressing the baby

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>



CARACTERISTICAS:

- Being a birthing baby he is flexible enough to fold into a pregnant belly or uterus and durable enough to withstand c-sections or vaginal birthing.
- The baby has lungs and can be intubated and ventilated with visible chest rise.
- His mouth and nose can be cleaned out.
- His umbilical stump has 3 bloodlines and blood can be taken from it. (Blood needs to be inserted before each use.)
- Each baby is unique in colouring and has eyelashes and curly black hair.
- Since the hair is very vulnerable it is wetted down with silicone to create the newborn wet look
- He can be washed, bathed, dried and powdered
- Length: 50cm. Head circumference: 37cm
- Weight: 3,2kg

INCLUI:

- 1 Baby Rohan
- 1 Set of clothes
- 1 Care kit: Baby powder, baby shampoo and dummy

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>