


# VASCULAR ACCESS ULTRASOUND PHANTOM - VAUS

Nº DE PRODUTO: 1891



## DESCRIÇÃO:

Vascular Access Ultrasound Phantom is an excellent healthcare simulator for the introduction and improvement of the techniques and psychomotor skills associated with successful ultrasound-guided vascular accessing.

The unique patient simulation material is flesh colored with the tactile feel and puncture resistance of human tissue.

## HABILIDADES:

- Ultrasound transducer positioning
- Vascular access skills training

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Espanha

### **CONTACTO**

sales@morethansimulators.com  
(+351) 21 123 85 18  
<https://morethansimulators.com>



## CARACTERISTICAS:

The curved contour on the top of the phantom is similar to contours on patients and has the added benefit of showing the importance of transducer positioning, as positioning can affect the quality of the image viewed—this is not possible on phantoms that are flat. This phantom can be imaged and accessed from the top, side or bottom surface, on the base or removed from the case.

The very durable ultrasound phantom can be accessed 1,000s of times, is self-healing and can be refilled with the supplied stain resistant simulated blood.

The Vascular Access Ultrasound Phantom has four vessels contained within the phantom:

- Very small superficial vein
- Mid-sized vein at medium depth
- Mid-sized vein which bifurcates into two smaller vessels
- Larger deeply place vein

## INCLUI:

- 1 Venous pad with 4 veins
- 1 Durable case
- 1 Simulated blood bottle

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Espanha

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+351) 21 123 85 18  
<https://morethansimulators.com>