

ECHOBABY - Newborn vascular access on arm

Nº DE PRODUTO: 2899



DESCRIÇÃO:

ECOBABY Vascular access arm phantom is an excellent healthcare simulator for the introduction and improvement of the techniques and psycho motor skills associated with successful ultrasound-guided vascular accessing in newborns and premature babies.

The unique patient simulation material is flesh colored with the tactile feel and puncture resistance of human tissue.

HABILIDADES:

- Ultrasound transducer positioning
- Newborn's vascular access skills training

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>

CARACTERÍSTICAS:

GENERAL:

- The curved contour of the phantom realistically reproduces the shape and anatomy of a newborn's arm.
- Is similar to contours on patients and has the added benefit of showing the importance of transducer positioning, as positioning can affect the quality of the image viewed–this is not possible on phantoms that are flat. This phantom can be imaged and accessed from the top, side or bottom surface, on the base or removed from the case.
- The very durable ultrasound phantom can be accessed 500s of times, is self-healing and can be refilled with the supplied simulated blood.

VEIN ACCESS:

- The ECOBABY phantom has two vessels contained within the phantom:
 - Small vein at medium depth
 - Very-small vein at medium depth

INCLUI:

- 1 ECOBABY task trainer
- 1 Syringe
- 1 Concentrated simulated blood
- 1 Rugged case PeliCase

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>