

SONOTRAIN - VENIPUNCTURE MODEL FOR ULTRASOUND

Nº DE PRODUTO: 2278



DESCRIÇÃO:

Ultrasound block with three blood vessels of diameter 4, 8 and 15 mm and adjustable flow rate.
The blood vessels are 30 to 40 mm below the surface.

HABILIDADES:

- Improve dexterity
- Improve hand-eye skill coordination
- Palpation and exploration skills
- Injection exercises

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>

- Venipuncture exercises

CARACTERISTICAS:

- Ultrasound block with three blood vessels of diameter 4, 8 and 15 mm and adjustable flow rate.
- The blood vessels are 30 to 40 mm below the surface.
- Possibility of puncture and injection, self-closing insertion channels.
- The SONOtrain™ trainer material simulates the sensation of real soft tissues to carry out palpations, injections, displays realistic texture and echogenicity in ultrasound images.
- It is an excellent teaching tool for basic training, to improve dexterity and hand-eye skill coordination.
- Ultrasound device sold separately
- Dimensions: 36 x 36 x 12,5 cm

INCLUI:

- 1 SONOtrain block
- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>