


WOUND PACKING AND BLEEDING CONTROL TRAINER

Nº DE PRODUTO: 2462



DESCRIÇÃO:

This Task Trainers is an excellent way to introduce your trainees and students to **realistic hemorrhage control for TCCC, TECC, TEMS, PHTLS**, and general wound packing techniques for advanced hemostatic gauze training and bleeding control techniques.

HABILIDADES:

- Bleeding point detection
- Wound packing skills training
- Hemostatic agents usage training
- Emergency and crisis management

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Espanha

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>



CARACTERÍSTICAS:

GENERAL:

- Upper realistic soft wound tissue and broken bone to add difficulty level.
- The task trainer is made from sturdy silicone, has no moving parts and is designed to be used over and over again.
- The integrated bleeding tubes can be connected to a manual or automatic pump.
- The leg can be washed and it is robust and reusable.

INCLUI:

- 1 Bleeding control task trainer with integrated tubes. Represents an open fracture.
- 1 Simulated blood
- 1 Filling syringe
- 1 Rugged case

MODULOS ADICIONAIS:

- TRUECLOT - Hemostatic gauze (sold separately)
- TRUECLOT - Cloggable simulated blood (sold separately)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+351) 21 123 85 18
<https://morethansimulators.com>