

# ADULT HEART & LUNG SOUNDS SKILLS TRAINER

Nº DE PRODUTO: 1221



## DESCRIÇÃO:

The HAL is a manikin designed to help **teach site-specific heart and lung auscultation skills** using seamless simulation technology. Simply press the bell against the skin and listen.

The students use realistic and palpable anatomic landmarks on the manikin to locate the site where a specific sound should be heard. If correct, they will hear the appropriate sound through the stethoscope and/or external speakers.

## HABILIDADES:

- Cardiac auscultation assessment
- Pulmonary auscultation assessment
- Stethoscope use

## CASOS E PATOLOGIAS:

### SOUNDS LIST:

- Base sound
- Physiological Split S2

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

- Paradoxical Split S2
- Fixed Split S2
- Split S2
- Opening snap
- Mid-Systolic Click
- Tracheal sounds
- Bronchial sounds
- Wheezing sounds
- Apex sounds
- Ronchi sounds
- Coarse crackles
- Med-Fine crackles
- Friction Rub
- Intermittent S4
- Stridor sounds
- Pleural friction
- Pulmonary Edema
- S3
- Starr-Edwards Valve

## CARACTERÍSTICAS:

### GENERAL:

- Adult size mannikin
- Hear the appropriate heart or lung sound as bell of stethoscope is moved across the front and back of the torso
- Full size adult torso with palpable anatomic landmarks
- Sensor network hidden beneath the skin
- Includes our Virtual Stethoscope® with multiple heart and lung sounds
- An external speaker can plug into the Virtual Stethoscope so a classroom can hear what the student hears (speaker not included)

## INCLUI:

- 1 HAL Skills trainer
- 1 Virtual stethoscope
- 1 Carrying bag
- 1 Instructions manual

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>