

# PRESTAN AED ULTRATRAINER

Nº DE PRODUTO: 2266



## DESCRIÇÃO:

The **Prestan AED UltraTrainer** offers features not found on any other universal AED trainer and is the most accurate representation of how current AEDs in the marketplace.

## CASOS E PATOLOGIAS:

- Includes **5 training scenarios**:
- Scenario 1 &gt; Shock / Shock / Shock(default)
- Scenario 2 &gt; Shock / No Shock / NoShock
- Scenario 3 &gt; No Shock / No Shock / NoShockd
- Scenario 4 &gt; Shock / No Shock / Shocke.Scenario 5 &gt; Shock / Shock / NoShock

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Espanha

### CONTACTO

sales@morethansimulators.com  
(+351) 21 123 85 18  
<https://morethansimulators.com>



- The AED trainer **can work as a Semi-automatic mode (SAED) or an Automatic Mode (AED):**
- Semi-Automatic &gt; Student must push flashing Shock Button to deliver shock(default)
- Automatic &gt; Unit will deliver shock when Shock Button is lit
  
- The trainer allows to configure a **time space to perform ventiliations**. This option can be activated or deactivated with a button.
- Ventilations On &gt; Time is provided in scenarios to give breaths(default)
- Ventilations Off &gt; No time is provided in scenarios to givebreaths

## CARACTERISTICAS:

### GENERAL:

- This compact and lightweight unit will give students a realistic training experience with training pads that sense when they are placed on a manikin.
- In addition, the pads are preconnected and voice prompts are clear and calm.
- Instructors can customize their training classes with options to turn on or off the compression metronome and/or the “give breaths” voice prompt.
- With a 3-year warranty and durable long-lasting pads, the AED UltraTrainer was built to last.
- Programming can toggle between two languages and is easily upgradeable for any AHA guideline changes. Comes standard with English/Spanish.
- 2015 AHA Guidelines compliant and updatable
- With play/pause button

### HIGH PERFORMANCE TRAINING PADS:

- Clear Anatomical Placement Graphics
- Long-Life Adhesive
- Adheres well to all Manikins with no messy residue left behind
- Unique Pad Connector with locking Connector Pin

### INTUITIVE AND COMPLETE INTERFACE:

- Prestan is the only universal AED Trainer with:
- Dual Graphic (Adult/Child) Training Pads
- Pad Sensing System

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Espanha

### CONTACTO

sales@morethansimulators.com  
(+351) 21 123 85 18  
<https://morethansimulators.com>



- Preconnected Pads
- Child Training Button
- Child AED specific voice prompts
- Automatic or Semi-Automatic shock simulation
- CPR Coaching, "Give Breaths" voice prompt
- Turn on or off compression metronome
- Turn on or off "Give Breaths" prompts
- Low Battery warning

## INCLUI:

- 1 AED trainer main unit
- 1 Pair of adult/paedriatic pad
- 1 Set of preinstalled batteries
- 1 Set of extra buttons
- 1 Carrying bag
- 1 Instructions manual

## MODULOS ADICIONAIS:

- **PAD REPLACEMENTS (sold separately):**
- 1 Pair of adult/child compatible pad replacement
- 4 Pairs of adult/child compatible pad replacements

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Espanha

### **CONTACTO**

sales@morethansimulators.com  
(+351) 21 123 85 18  
<https://morethansimulators.com>