

# AirSim Baby X

Nº DE PRODUCTO: 1101



## DESCRIPCION:

TruCorp AirSim Baby X es el entrenador de vías respiratorias pediátrico anatómicamente más correcto para enseñar y capacitar en el manejo de las vías respiratorias infantiles. Los detalles anatómicos internos y externos precisos, basados en datos de CT-DICOM de un niño de seis meses, proporcionan una práctica realista para las actividades pediátricas de intubación, laringoscopia, ventilación e inserción de dispositivos supraglóticos.

## HABILIDADES:

- Double nasotracheal intubation
- Bag and Mask ventilation techniques
- Full range of supraglottic devices
- Direct and video laryngoscopy

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05  
<https://morethansimulators.com>

- Endotracheal tube insertion

## CARACTERÍSTICAS:

- AirSim X airway is certified to practice 20,000+ intubation cycles without fail and is supported by a 5-year airway guarantee
- User feedback through the visualization of stomach inflation and lung expansion
- Modular design allowing for all parts to be changed within your organization (if required)
- Delivered ready to use in a durable carry case for easy transportation and safe storage

## INCLUYE:

- 1 AirSim X (ensamblado y listo para usar)
- 1 Trucorp Carry Case
- 1 USB PEN Drive que contiene manual del usuario del producto
- 1 botella de lubricación de Trucorp (100 ml)

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05  
<https://morethansimulators.com>