


# PEDI Cuidados pediátricos avanzados un año

Nº DE PRODUCTO: 1283



## DESCRIPCION:

One-Year-Old PEDI® Nursing Care Patient Simulator permite a los participantes practicar y desarrollar habilidades de cuidados de enfermería neonatal. Incluye una vía aérea intubable, un brazo de entrenamiento de inyección y una pierna de inyección intraósea para que los participantes puedan trabajar de manera práctica durante los ejercicios de compresión torácica y ventilación.

## HABILIDADES:

- Bandaging
- Hygiene care
- Palpation of pulses
- Ear canal examination
- BVM ventilation

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05  
<https://morethansimulators.com>

- Airway management
- Nasal intubation
- Intubation, suctioning and Sellick maneuver
- NG and OG tube feeding
- Gastric suction
- Intramuscular injection exercises
- Intraosseous infusion
- IV training arm with intramuscular, intradermal and subcutáneo

## CARACTERÍSTICAS:

### FEATURES:

- Full-body with soft, lifelike face skin and molded hair
- Eyes can be closed or opened manually
- Head tilt/chin lift
- Endotracheal intubation
- Nasogastric intubation
- Tongue and teeth
- Intraosseous infusion
- IV training arm
- Femoral vein
- Manual palpable pulses: femoral and brachial
- BVM with realistic chest rise
- CPR anatomical landmarks
- Lavage/gavage
- Intramuscular injection sites
- Simulated ear canal
- Practice gastric suctioning exercises
- Articulating head, arms, and legs
- Realistic hands, feet, fingers, toes

## INCLUYE:

- 1 maniquí pediátrico de tamaño completo
- 1 conjunto de ropa (camiseta + pantalones cortos)
- 1 bolsa de transporte suave
- 1 Manual de instrucciones

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05  
<https://morethansimulators.com>