

# TENSION PNEUMOTHORAX SIMULATOR

Nº DE PRODUCTO: 2298



## DESCRIPCION:

Recommended by EMT (Emergency Medical Training) instructors, this model permits instruction of the delicate procedure for **successfully managing chest wounds in which a collapsed lung interferes with the victim's respiration and blood flow.**

## HABILIDADES:

- Needle insertion in second or fifth intercostal space
- Thoracotomy in second or fifth intercostal space

## CARACTERISTICAS:

- This model allows medical personnel to practice proper needle insertion which enables built-up air pressure to escape the thorax.
- The model accepts needle and thoracotomy in second intercostal space in the midclavicular line or fifth intercostal space in the midaxillary line.
- Dimensions: 23 x 48 x 48 cm
- Weight: 3.5 kg

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05  
<https://morethansimulators.com>



**More Than Simulators™**  
*Excellence. Every patient. Every time.*

## INCLUYE:

- Tension pneumothorax simulator
- 1 Rib cage
- 1 Foot pump
- 1 Pleural chamber replacement
- 1 Soft carrying case
- Instructions manual



### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05  
<https://morethansimulators.com>