

EAR MODEL FOR TYMPANOPLASTY

Nº DE PRODUCTO: 2114



DESCRIPCION:

Este modelo de oreja de silicona está diseñado para ayudar a enfermeros, asistentes y profesionales de la salud a entrenar varios ejercicios de entrenamiento auditivo y procedimientos quirúrgicos como la extracción de objetos extraños y la timpanoplastia.

HABILIDADES:

- Ear canal foreign objects removal
- Anaesthesia administration
- Tympanoplasty
- Ear lobe piercing
- Ear canal irrigation techniques

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05
<https://morethansimulators.com>

CARACTERÍSTICAS:

GENERAL:

- The model is made from high-quality silicone.
- The base of the model has a rugged surface to prevent the model from gliding.
- The eardrum membrane is easily replaceable to allow multiple tympanoplasty training sessions.
- The lobes can be perforated to train cosmetic and surgical exercises

INCLUYE:

- 1 modelo de entrenamiento de oído
- 20 membranas de tímpano

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05
<https://morethansimulators.com>