



# ORAL ANESTHESIA MANIKIN

Nº DE PRODUCTO: 3041



## DESCRIPCION:

Practice local anesthesia techniques.

## HABILIDADES:

- Dry injection
- Standard Mandibular block anesthesia
- Gow Gates anesthesia
- Local infiltrations

## CARACTERISTICAS:

### FEATURES:

- Model consists of a metal skull with flexible silicone face and mouth lining over an anatomically correct jaw with 28 permanent teeth, movable tongue, and jawbone with landmarks for placement of dry injection that can be visibly seen and palpated.
- Conduction sensors aid in practicing Gow Gates and Standard Mandibular block anesthesia.
- Local infiltrations may also be administered in their respective regions
- Battery powered LED light on top of the manikin's head indicates when the correct position and angle of injection is achieved by signaling

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05  
<https://morethansimulators.com>



**More Than Simulators™**  
*Excellence. Every patient. Every time.*

a green light.

- Weight: approx 4kg

## INCLUYE:

- 1 Oral anesthesia manikin

- 3 Batteries

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05  
<https://morethansimulators.com>