

# PAT - PEDIATRIC AUSCULTATION TRAINER

Nº DE PRODUCTO: 1388



## DESCRIPCION:

PAT – the Pediatric Auscultation Trainer – is the complete auscultation trainer for pediatric-specific conditions.

Because pediatric care requires careful consideration of growth and developmental issues, PAT is designed to prepare students to readily recognize the signs and sounds associated with pediatric-specific disorders.

## HABILIDADES:

- Auscultation of hear, lung, bowel and bruit sounds
- Stethoscope usage skills
- Phonocardiograms monitoring and diagnosis

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05  
<https://morethansimulators.com>



## CARACTERÍSTICAS:

### KEY FEATURES:

- With one of the largest libraries of pediatric sounds and intuitive instructional tools, PAT gives students everything they need to actively engage in their auscultation education.
- Large library of pediatric-specific heart, lung, bowel and bruit sounds
- Students use their own stethoscope (clean and safe)
- Variable breath sounds between left/right lobes
- Phonocardiograms
- Smart classroom/auditorium ready
- Customized volumen adjustment
- Lessons guide for each condition
- Customizable lectures and scenarios

### PEDIATRIC-SPECIFIC SOUNDS:

- Effective assessment of pediatric patients is critical for ensuring proper care of this rapidly developing population.
- With that in mind, we built a library of sounds tailored to pediatric curricula – whether secondary or graduate level.
- Now includes a COVID-19 sound.

### CLEAR, ACCURATE SOUNDS:

- Providing one of the largest libraries of pediatric sounds, PAT produces the clearest, most accurate sounds at their anatomically-correct locations.
- Students can auscultate anteriorly and posteriorly from a collection of over 100 cardiac, respiratory, bowel, and bruit sounds.
- And with built-in speakers, students are able to use their own stethoscopes – helping build confidence in their skills.

### ULTRA PORTABLE:

- PAT is incredibly lightweight and portable.
- Easily integrated into any smart classroom or auditorium, PAT can enhance lectures by projecting sounds aloud and by modeling proper auscultation methods.

### COMPLIMENTARY 1-YEAR SAM ONLINE LICENSE:

- PAT users receive a complimentary license for one year to SAM Online – an online platform with an extensive sounds library, comprehensive lesson guides, and tools to examine progress around-the-clock.
- Together with PAT, students will have constant access to the tools they need to master their auscultation skills.

## INCLUYE:

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05  
<https://morethansimulators.com>

- 1 PAT Auscultation manikin
- 1 Laptop with preinstalled software
- 1 Digital copy Lessons guide (English, Spanish or German)
- 1 Operators manual

## MODULOS ADICIONALES:

### **SAM ONLINE (1 year access included):**

- AM Online is an online platform that gives students and instructors access to an extensive simulated sounds library, lesson guides for over 60 conditions, and tools to examine progress around-the-clock.
- Available for one year to users of SAM, PAT, and the Bionic Hybrid Simulator – or for individual purchase – SAM Online is an essential program for reinforcing the concepts of auscultation at any academic level.
- Now includes COVID-19 sound.
- A full license includes 1 admin license, 5 teacher licenses, and 200 student licenses.

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05  
<https://morethansimulators.com>