

SIMULADOR DE PATA CANINA PARA ACCESO VASCULAR

Nº DE PRODUCTO: 2204



DESCRIPCION:

Entrene técnicas de venopunción y acceso vascular con nuestro Canine Leg Vascular Access Simulator.

CARACTERISTICAS:

- Made of high quality durable silicones, the Canine Leg is designed specifically for blood collection and catheter insertion of a 22 gauge or smaller.
- Our simulated cephalic vein is softer than typical latex tubing to enhance realism and effectiveness.
- This model gives students the ability to practice their vascular access techniques without the need for a partner or a live animal.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05
<https://morethansimulators.com>



More Than Simulators™
Excellence. Every patient. Every time.

INCLUYE:

- 1 Simulador de acceso vascular de la pierna canina
- 1 bolsa IV
- 1 Polvo de sangre artificial
- 1 reemplazo de la piel
- 1 Manual de instrucciones



MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05
<https://morethansimulators.com>