

# AirSim Child Bronchi X

Nº DE PRODUCTO: 1104



## DESCRIPCION:

TruCorp AirSim Child Bronchi X ofrece a médicos pediátricos, especialistas clínicos pediátricos y equipos de emergencia una forma de practicar competencias de manejo de las vías respiratorias pediátricas y técnicas de diagnóstico en un maniquí anatómicamente correcto. Adecuado para la formación en laringoscopia, broncoscopia, intubación y técnicas de aislamiento pulmonar. Incluye la vía aérea AirSim X con una garantía de cinco años.

## HABILIDADES:

- Oral and nasal intubation
- Bag/mask ventilation
- Laryngeal Mask Airway (LMA) and supraglottic device insertion
- Combitube insertion; Double lumen endotracheal tube insertion

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05  
<https://morethansimulators.com>

- Diagnostic bronchoscopy
- Lung isolation techniques using left and right endobronchial tubes and bronchial blockers
- Lung suctioning techniques

## CARACTERÍSTICAS:

- AirSim X airway is certified to practice 20,000+ intubation cycles without fail and is supported by a 5-year airway guarantee
- Features bronchi which provides internal anatomically correct features down to the fourth generation
- Inflatable tongue bulb to create tongue edema providing a variance during training
- User feedback through the visualization of stomach inflation and lung expansion
- Missing front tooth to illustrate realistic dentals of a 6-year old
- Modular design allowing for all parts to be changed within your organization (if required)
- Delivered ready to use in a durable carry case for easy transportation and safe storage

## INCLUYE:

- 1 AirSim Child Bronchi X (ensamblado y listo para usar)
- 1 Trucorp Carry Case
- 1 USB PEN Drive que contiene manual del usuario del producto
- 1 botella de lubricación de Trucorp (100 ml)

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05  
<https://morethansimulators.com>