

SUPER OB SUSIE - Postpartum hemorrhage kit

Nº DE PRODUCTO: 1670



DESCRIPCION:

Postpartum hemorrhage package enables the training of postpartum emergency management skills.

HABILIDADES:

- Practice catheterization with the urethra that connects to a realistic urinary bladder.
- Adjust uterine tone with the pump attachment
- Perform fundal massage.
- Practice Bakri balloon tamponade placement.
- Rectal suppository administration.

CARACTERISTICAS:

FEATURES:

Postpartum hemorrhage capability is available through the blood reservoir that holds up to 400mL of fluid.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05
<https://morethansimulators.com>



INCLUYE:

- 1 Postpartum hemorraghe perineum insert
- 1 Boggy uterus
- 1 Postpartum hemorraghe cervix
- 1 Elevating pillow
- 1 Blood concentrate
- 1 Fill kit
- 1 Instructions manual



MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05
<https://morethansimulators.com>