

# HAL HEART & LUNG SOUNDS SKILLS TRAINER

PRODUCT NUMBER: 1221



## DESCRIPTION:

The HAL® is a manikin designed to help teach site-specific heart and lung auscultation skills using seamless simulation technology. Simply press the bell against the skin and listen.

The students use realistic and palpable anatomic landmarks on the manikin to locate the site where a specific sound should be heard. If correct, they will hear the appropriate sound through the stethoscope and/or external speakers.

## SKILLS:

- Cardiac auscultation assessment
- Pulmonary auscultation assessment
- Stethoscope use

## CASES AND PATHOLOGIES:

### SOUNDS LIST:

- Venous hum
- Aortic stenosis

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- Split S2
- Pulmonary stenosis
- Systolic fixed S2
- 1 year heart
- 6 year heart
- Stills murmur
- Mitral valve regurgitation
- Split S1
- Normal infant
- Normal child
- Stridor sounds
- Wheezing sounds
- Ronchi sounds

## CHARACTERISTICS:

### FEATURES:

- Adult size mannikin
- Hear the appropriate heart or lung sound as bell of stethoscope is moved across the front and back of the torso
- Full size adult torso with palpable anatomic landmarks
- Sensor network hidden beneath the skin
- Includes our Virtual Stethoscope® with multiple heart and lung sounds
- An external speaker can plug into the Virtual Stethoscope so a classroom can hear what the student hears (speaker not included)

## INCLUDE:

- 1 HAL Skills trainer
- 1 Virtual stethoscope
- 1 Carrying bag
- 1 Instructions manual

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>