

HAL AUSCULTATION, AIRWAY & CPR SKILLS TRAINER - OMNI2

PRODUCT NUMBER: 1220



DESCRIPTION

HAL® is an easy-to-use portable airway management and CPR skills trainer designed for teaching resuscitation and auscultation skills to students and professionals in EMS, Medicine, and Nursing. HAL® features a high-fidelity airway, realistic heart and lungs sounds, and real-time CPR quality feedback to help you teach effective cardiopulmonary resuscitation and airway management.

HAL® is also available in a full-body configuration.

SKILLS

- Airway management
- Nasal and oral intubation
- Resuscitation training (CPR)
- BVM ventilation
- Cardiac and pulmonary auscultation assessment

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 671 222 375 / (+351) 917 260 661
www.morethansimulators.com



Debriefing tool

CHARACTERISTICS

FEATURES:

Adult-sized, male upper torso

Anatomically correct landmarks for proper CPR hand placement

AIRWAY MANAGEMENT:

Realistic airway with tongue, vocal cords, and trachea

Realistic head tilt/chin lift and jaw thrust

Nasal and oral Intubation with standard adjuncts: ETT, LMA, King LT®

CPR hand position landmarks

Visible chest rise with bag-valve-mask ventilation

Esophageal intubation and/or excessive ventilation produce gastric distension

Right mainstem intubation presents unilateral chest rise

CPR:

Compliant with 2015 adult CPR guidelines

Realistic chest cavity resistance and recoil

Real-time feedback of CPR quality: Compression depth and rate, Ventilation Rate, Excessive Ventilation, No-Flow Time and CPR Cycles

CARDIAC/CIRCULATION:

Palpable carotid pulses ? Manually generated

40 Programmable abnormal and normal heart and lung sounds

INCLUDE

1 HAL Skills trainer

1 OMNI2 Wireless Control Interface (Tablet+License)

1 Carrying case

1 Power supply adapter

1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 671 222 375 / (+351) 917 260 661
www.morethansimulators.com