




VASCULAR MICRO-SURGERY HEAD SIMULATOR

PRODUCT NUMBER: 2449



DESCRIPTION:

The **Carotid Head Simulator** consists of a male size head with a cervical vascular wound at the neck where a Carotid Artery bifurcation can be placed to provide a highly realistic vascular micro-surgery training environment.

SKILLS:

- Carotid vascular surgery techniques
- Eye-to-hand coordination
- Surgical instruments handling and recognition
- Carotid trauma open repair procedure
- Carotid endarterectomy Artery (CAE) procedure
- End-to end and end-to-side anastomosis training

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



CHARACTERISTICS:

GENERAL:

- Size: 290 x 200 x 90 mm
- Weight: 350 gr

FEATURES:

- This ultra-realistic training simulator is the perfect task trainer for those who want to test and improve their surgical skills abilities in a true-to-life environment. From novel surgeons to experienced professionals, this simulator can adapt to all levels of medical proficiency.
- Designed to be used with real instrumental and equipment.
- Includes integrated bleeding tubes. When connected and filled to a pump, can recreate the human bleeding.
- This simulator can be customized by adding simulated plaque to the arteries to increase the level of the training.
- Manufactured with high-grade silicone and latex. It avoids the necessity to acquire animal tissue and its difficult storage conditions.

INCLUDE:

- 1 Carotid Head Simulator
- 1 Stand
- 1 Set of carotid arteries
- 1 Simulated Blood

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>