




# ZACK - Multipurpose male care

PRODUCT NUMBER: 1198



## DESCRIPTION

The ZACK? Multipurpose Male Care Simulator is a **full-sized, adult male lower torso** (abdomen and pelvis) designed to provide health care professionals with hands-on **training in the skills required to perform men's healthcare procedures**.

It is useful for demonstrating and practicing catheterization exercises, prostate examinations, testicular examinations, no-scalpel vasectomy, and rectal/colon examinations.

## SKILLS

- Catheterisation exercises
- Prostate examination
- Testicular examination
- No scalpel vasectomy
- Rectum and colon examinations



## CASES AND PATHOLOGIES

Melanoma  
Ulcerative adenocarcinoma  
Adenomatous polyp  
Sessile adenocarcinoma  
Advanced polypoid adenoma  
Pedunculated adenomatous polyp  
Pedunculated polyp

## CHARACTERISTICS

Full-size male lower torso  
Internal bladder for catheterization  
Four interchangeable prostates depicting a moderately enlarged benign prostate, a prostate with two discreet nodules, a prostate with an easily palpable large mass, and a prostate with malignant cancer  
Realistic penis and two scrotal sacks. One scrotal sack is normal, and the other contains tumors in each testicle  
No-Scalpel Vasectomy kit containing two removable scrotal skins, two testicles, and two long vas assemblies that can be advanced as needed for NSV exercises  
Rectum and colon containing benign and malignant masses easily visualized using an appropriate endoscope

## INCLUDE

1 Lower torso male simulator  
1 Internal bladder for catheterization  
4 Interchangeable prostates  
1 No scalpel vasectomy kit  
1 Rectum and colon with benign and malignant masses  
Talcum powder  
1 Soft carrying bag  
1 Instructions manual