



SkinPad - Incision and suture skills trainer

PRODUCT NUMBER: 1017



DESCRIPTION

Completely new **3 layer skin pad** for demonstrating and practising a wide variety of incisions and suturing skins.

It is available in **medium size** (120x75mm) or **large size** (120x150mm).

SKILLS

- Incisions: linear, ellipse, flaps and shaped
- Intradermal suturing
- Simple and advanced interrupted suturing
- Subcuticular suturing
- Continuous suturing
- Stapling



Use of adhesive strips

CHARACTERISTICS

Dimensions: 120 x 75 mm o 120 x 150 mm

Realistic representation of the 3 main skin layers: epidermis, dermis and subdermal layer (fascia).

Realistic retention of sutures thanks to its high-quality silicone

Integrated mesh to increase the resistance of the tissue

Adequate for beginners and specialists

Reusable. The entire surface allows incisions and anaesthetic injection.

It is recommended to use for one year after purchase. High quality silicone includes oils in its manufacture that guarantee a correct texture. These oils evaporate over time.

INCLUDE

1 SkinPad Thin (can be ordered in packs of 2, 4, 6 or 10 units)

ADDITIONAL MODULES

INTRUMENTS CASE (option): Add a 19-pieces dissection kit to enhance your suture practice.

LIPOMA PAD (sold separately): A realistic pad with 2 lipoma to practice full lipoma removal.

SEBACEOUS CYST PAD (sold separately): A realistic pad with 3 cysts to practice full sebaceous cyst removal procedure.