

REALITi Plus - Simulated patient monitor

PRODUCT NUMBER: 1885



DESCRIPTION:

The iSimulate REALITi Plus Simulated Patient Monitor is the basic version of the REALITi360 ecosystem and enables medical teachers to reproduce multiple scenarios, whether basic or complex.

It is a true mobile simulation center housed in a bag that gives you the freedom to work anywhere, because it allows you to train and evaluate directly in ambulances, helicopters, hospital settings, simulation centers or training laboratories.

CASES AND PATHOLOGIES:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



STANDARD SCREENS:

- Defibrillator screen
- Monitor screen
- DEA screen
- ICU 1 screen
- ICU 2 screen

ADVANCED SCREENS:

- Corpuls1 patient monitor screen
- Corpuls3 patient monitor screen
- Corpuls® AED/defibrillator screen
- CARESCAPE™ B40 patient monitor screen
- LIFEPAK® 1000 patient monitor screen
- LIFEPAK® 15 patient monitor screen
- LIFEPAK® 20 patient monitor screen
- Medtronic Capnostream™ 35 patient monitor screen
- Philips IntelliVue MX800 patient monitor screen
- Philips HeartStart MRx emergency care patient monitor screen
- Philips HeartStart MRx hospital patient monitor screen
- Welch Allyn Connex® VSM 6000 patient monitor screen
- Zoll® Propaq® MD patient monitor screen
- Zoll® R Series® patient monitor screen
- Zoll® X Series® patient monitor screen
- Schiller DEFIGARD Touch 7 patient monitor screen
- Schiller PHYSIOGARD Touch 7 patient monitor screen

CHARACTERISTICS:

FEATURES:

Thorough emphasis on detail.

Same look, sound, and performance as a real patient monitor.

Separate respirator simulator that can be used separately.

The respirator works in both portrait and landscape orientation.

VERSATILITY:

The user can add ectopic pregnancies, arrhythmias, and artifacts to any ECG wave.

Any ECG wave can have premature ventricular contractions (CVP), premature atrial contractions (CAP), and premature junction complexes (CUP).

Sinus arrhythmia linked to respiratory rate.

Artifacts and electrical interference.

Visualization of trends in the wave.

LAB AND RADIOLOGY DATA-BASE:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



Radiological images processed

CT and X-rays

Laboratory results with more than 100 parameters

LEARNING MANAGEMENT SYSTEM:

Event log, CPR compressions, vital signs history, score, waves, video and all other data stored in Control.

The video is transferred to Control at the end of a scenario and can be exported with overlapping vital sign parameters.

Exploring students and previous scenarios.

Export of scenario reports in PDF format for sharing.

Export of all data in CSV format for sharing.

INCLUDE:

- 1 Monitor ipad with pre-installed REALITi360 app
- 1 Facilitator ipad with pre-installed REALITi360 app
- 1 Set of accessories:
 - ECG cable (3 derivations)
 - ECG cable (7 derivations)
 - Troncal ECG cable (10 derivations)
 - SpO2 sensor
 - Adult PNI cuff
 - DEA/Defibrillator connector
 - Disposable ECG adult electrodes
- 1 Wifi Router
- 1 Carrying bag
- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>